



emulate

Protocol for Emulate Organ-Chips:
Cytokine-Mediated Barrier Disruption

April 8, 2021

EP213 v1.0

TITLE Cytokine-Mediated Barrier Disruption	DOCUMENT EP213	VERSION 1.0
	DATE 4-8-2021	PAGE 2 OF 4

Goals:	Key Steps:	Other Required Materials:
<ul style="list-style-type: none"> • Prepare treatment solution • Treat cells in Emulate Organ-Chips with IL-22 	<ul style="list-style-type: none"> • Prepare stock solution • Prepare working solution • Treat chips with cytokine 	<ul style="list-style-type: none"> • IL-22 (R&D, Cat. # 782-IL-010) • IL-22BP(R&D, Cat. # 8498-BP-025) • IFNγ (Peprotech, Cat. # 300-02) • Tofacitinib (Sigma, Cat. # PZ0017) • DMSO (Sigma, Cat. # D2650) • BSA (Sigma, Cat. # A1595) • DPBS (Sigma, Cat. # D8537) • Medium

Introduction

The Colon Intestine-Chip provides a human organoid-based platform to study cell-cell interactions and mechanisms that elicit collapse of the intestinal epithelial homeostasis. We employed the cytokines IFN γ , a prototype cytokine on intestinal epithelial barrier disruption studies, or IL-22. IL-22 acts as a barrier-insulting cytokine when applied to healthy tissue in our Colon Intestine-Chip, in line with recent findings^{1,2,3}. For barrier disruption of customer-originated cell material, this protocol may need to be optimized. Please reach out to [Emulate Field Science Support](#) for additional guidance.

Method

Before starting this process, ensure that all treatment solutions prepared in Part I and II of this protocol are warmed up to 37°C. Label all chips and Pods with their corresponding treatment conditions prior to adding the compound.

Part I - Barrier function mediated by IL-22

1. Prepare the stock solution, diluting the cytokine in the appropriate solvent (e.g., DMSO, water, DPBS), as per manufacturer's instructions. Ensure that all materials and reagents required for the preparation of the stock solution are sterile and ready.
2. Calculate the volume of solvent to prepare a stock solution of the test article based on mass / molecular weight / molarity of the compound and the desired concentration / molarity of the stock solution. Preparation of 1000X target concentration of stock solution is typically recommended as this allows you to maintain solvent at 0.1%.
3. Resuspend IL-22 in DPBS at stock concentration of 10 μ g/mL (606.06 nM).
4. Resuspend the IL-22 BP in 0.1% BSA in DPBS at concentration of 250 μ g/mL (9,615.38 nM).
5. Ensure all chip carriers are labeled and identify the different conditions clearly.

TITLE Cytokine-Mediated Barrier Disruption	DOCUMENT EP213	VERSION 1.0
	DATE 4-8-2021	PAGE 3 OF 4

Part II — Barrier Function mediated by IFN γ

1. Reconstitute the IFN γ in sterile cell culture grade water. We recommend to resuspend the 20 μ g vial in 200 μ L of prepare the stock solution at a final concentration of 10 μ g/mL.
2. Once reconstituted, aliquot and store the stock solution in -80°C. Do not keep the aliquots for longer than one week.
3. Reconstitute 5 mg of Tofacitinib in 495.55 μ L of DMSO to prepare a stock solution at a concentration of 20 mM.
4. Aliquot the reconstituted Tofacitinib in 20 μ L portions. This amount is sufficient to dose 12 chips.

Part III – Dosing the Cells in the Colon Intestine-Chip

1. Carefully remove trays with Pods from the Zoë™ Culture Module and transfer to the biosafety cabinet. We recommend removing one tray at a time to minimize stress experienced by cells because of fluctuations in temperature. Ensure that the biosafety cabinet light is turned off if you are working with a light sensitive compound.

Note: If you have multiple collection time points in your experiment, it helps to organize your conditions so that you have one time point per tray. This will allow you to pause flow only on the tray from which you need to collect, while letting other trays flow uninterrupted, resulting in a more accurate assessment of elapsed time.
2. Fully aspirate both inlet and outlet reservoirs of each Pod, while avoiding direct contact with the Pod reservoir vias.
3. Add the calculated volume of warm, freshly-prepared treatment media to the appropriate channel.
4. Once all Pods have been refreshed, ensure that all trays are returned to the appropriate Zoë. Flush chips at 600 μ L / hour for 5 minutes to introduce dosing solution in the chip channels.
5. After flushing Pods and chips, carefully remove trays from Zoë one at a time. Once again, transfer to the biosafety cabinet, and aspirate from the flow-through medium in the outlet reservoirs.
6. Return all trays to Zoë, reset the Zoë settings to the correct experimental conditions (e.g., flow rate, stretch) and note the time as the experimental start time (t=0).
7. Each outlet reservoir can be sampled independently at each timepoint following Protocol EP124 Effluent Sampling.
8. Replenish Pods with freshly prepared treatment medium daily, regardless of the collection timepoints, until the end of your treatment period or experiment.

Test Article	Concentration	Starting Day	Treatment Duration	Flow rate
IL-22	1 nM	Day 5	1-72 hours	30 μ L/h
IL-22BP	30 nM	Day 5	1-72 hours	30 μ L/h
IFN γ	10 ng/mL	Day 4	48 hours	30 μ L/h
Tofacitinib	20 μ M	Day 4	72 hours	30 μ L/h

TITLE Cytokine-Mediated Barrier Disruption	DOCUMENT EP213	VERSION 1.0
	DATE 4-8-2021	PAGE 4 OF 4

Suggested Readouts:

- [RNA Isolation and Purification Protocol](#)
- [cDNA Preparation Protocol](#)
- [Quantitative Polymerase Chain Reaction \(qPCR\) Protocol](#)
- [Protein Samples Isolation Protocol](#)
- [Total Protein Quantification Protocol](#)
- [Assessment of Cleaved Caspase 3 Protocol](#)
- [Assessment of Phosphorylated STAT3 Protocol](#)
- [Quantification of Cytokines and Acute Inflammatory Phase Proteins Secretion Protocol](#)
- [Immunofluorescence Staining](#)

References

1. Wang, Y., Mumm, J. B., Herbst, R., Kolbeck, R. & Wang, Y. IL-22 Increases Permeability of Intestinal Epithelial Tight Junctions by Enhancing Claudin-2 Expression. *J. Immunol.* (2017). doi:10.4049/jimmunol.1700152
2. Powell, N. *et al.* Interleukin-22 orchestrates a pathological endoplasmic reticulum stress response transcriptional programme in colonic epithelial cells. *Gut* (2019). doi:10.1136/gutjnl-2019-318483
3. Tsai, P. Y. *et al.* IL-22 Upregulates Epithelial Claudin-2 to Drive Diarrhea and Enteric Pathogen Clearance. *Cell Host Microbe* (2017). doi:10.1016/j.chom.2017.05.009

© Emulate, Inc., 2019. All rights reserved.

Emulate® and the Emulate logo are registered trademarks of Emulate, Inc.

The technology disclosed in this document may be covered by one or more patents or patent applications owned by or licensed to Emulate, Inc. No license is granted herein. You are solely responsible for determining whether you have all intellectual property rights that are necessary for your intended use of Emulate products or protocols, and whether you are required to obtain any additional intellectual property rights from a third party. Further information is available by contacting Emulate.